EN

Easy preloaded IOL – Hydrophilic Acrylic, Preloaded Intraocular Lens

Description:

The Easy IOL is a preloaded intraocular lens supplied with a compatible injector. The lens is made of hydrophilic acrylic material (25% water content) with ultra violet blocker and violet filtering chromophore for increased protection of the retina.

The Easy IOL is supplied in one of the two platforms - BunnyLens (4-loop) or SeeLens (C-loop).

See label on the cardboard box for type of the lens:

Label	Optical design
AF	Monofocal Aspheric
MF	Multifocal Diffractive Apodized Aspheric <ul style="list-style-type: none">The label indicates the distance power. Addition for near vision is +3 diopters.
TR	Toric Aspheric. <ul style="list-style-type: none">The label indicates the Spherical Equivalent and Cylinder powers.

Indication:

The Easy IOL is indicated for surgical treatment of senile cataract, intended for placement in the capsular bag.

Contraindications:
Absolute contraindications: Any chronic condition where an undesirable outcome is expected. This may include:

- Chronic acute uveitis.
- Retinal diseases in which the implant may interfere with retinal surgery.
- Rubella cataract
- Progressive diseases of the anterior segment.

Relative contraindications:

- Significantly irregular corneal topography or previous corneal transplant.
- High astigmatism or keratoconus
- Severe corneal dystrophy
- Amblyopia
- Uncontrolled glaucoma
- Optic nerve atrophy
- Any case requiring intraoperative manipulation to enlarge the pupil.
- Aniridia or iris neovascularization
- Microphthalmals or macrophthalmos
- Intraoperative Hyphema, significant vitreous loss or bleeding
- Uncontrollable intraoperative intraocular pressure.

Clinical cases which may deteriorate due to IOL implantation and cases with an increased risk of implantation based on the surgeon's experience. The evaluation of each individual case is under the surgeon discretion.

Complications:
Cataract surgery, with or without lens implantation, might be associated with:

- Ocular inflammation
- Hemorragia
- Intraocular pressure elevation
- Post operative infection
- Retinal detachment
- Macular edema
- Corneal edema
- Posterior capsule opacification
- Capsular rupture
- Vitreous loss

Complications related to intraocular lens implantation:

- Lens decentration and luxation
- Inaccurate lens power calculation
- Damage to lens during implantation

Warnings:

- The IOL must be implanted in accordance to the following instructions for use. Improper use may impose risk to the patient's health.**
- Do not use the intraocular lens if the external sterile packaging is damaged, or when leakage from the container has occurred, or in any case of damage to the IOL. It is essential to wet it in balanced solution or equivalent, before implantation or when handing over a longer period of time.
- The lens might have a violet reflecting tint while inspecting through the slit lamp due to its transmittance properties, depending on slit lamp intensity and angle.
- Patients intended for refractive lens exchange may be at greater risk for retinal detachment. It is recommended to determine the relative benefit to the patient in these cases. Lens removal techniques such as low power pulsed ultrasonic phacemulsification and liquefaction may lower this risk, whereas the surgeon should avoid high power phacemulsification and Perts vitreus.
- The surgeon should target emmetropia to achieve optimal results.
- The lens should be implanted so that it is centered optimally, in order to achieve optimal results and to avoid visual disturbances.
- Patients with preoperative or expected postoperative astigmatism of >1.0D may not achieve an optimal visual outcome and satisfaction.
- A high level of surgical skill is required for implantation of intraocular lenses. It is recommended that the surgeon observe and assist in several procedures prior to attempting implantation.
- Special considerations for toric (TR) IOLs:
 - The lens should be implanted so that it is centered optimally, in order to achieve optimal results and to avoid visual disturbances.
 - The lens should be implanted following the instructions described in Hanita Lenses Toric IOL calculator printout.
 - In any case of capsular rupture, zonular damage or a posterior capsulotomy is planned, the lens should not be implanted.
 - It is important to remove any remaining viscoelastic material at the end of surgery, as any residues may cause unwanted rotation of the IOL leading to residual postoperative astigmatism.
 - A high level of surgical skill is required for implantation of intraocular lenses. It is recommended that the surgeon observe and assist in several procedures prior to attempting implantation.

Packaging:
The Easy IOL is supplied sterile in a blister pack inside a peel pouch. The Easy IOL is placed in a loading chamber, fixed by a lens holder, all immersed in saline.

The packaging sterility is guaranteed unless the bag is opened or damaged.
Instructions for use:

There are various surgical procedures which can be utilized. The surgeon should select a procedure which is appropriate for the patient.

- Examine the label on the lens package for IOL type, diopter and expiration date. In a sterile environment, open the peel pouch and remove the blister pack. Verify the diopter of the lens again.
- Inspect the blister pack. Make sure it is not damaged and the seal is not broken. Before opening the blister pack, gently tap on the lid to remove drops of storage liquid from the inside of the lid.
- Slowly and continually peel off the lid while holding the blister pack in a horizontal position.
- Remove the loading chamber containing the lens from the blister pack. Do not remove the lens holder at this time!
- Insert the loading chamber into the compatible injector. When inserting the loading chamber into the injector ensure that the plunger is in withdrawn position.
- Apply an appreciate amount of viscoelastic to the loaded injector. Ensure that the visco elastic migrates under the lens and fills the injector tip.
- Remove the lens holder from the loading chamber.
- Close the wings of the loading chamber. The system is ready for injection.

For detailed descriptions, refer to the instructions for use supplied with the compatible injector.

- Once the IOL is in place, thoroughly remove excess viscoelastic from the eye, both in front and behind the IOL, by routine irrigation and aspiration. Alternatively, the lens can be carefully removed from the loading chamber, taking care not to touch the optic part. The lens can be then be inserted using either forceps or conventional injector.

Liability:

Hanita Lenses company covers the design and the production of the intraocular lens. It shall be incurred in no way in case of accidents resulting from the use of this lens. The lenses are carefully checked and inspected by the manufacturer to assure a high quality product. If a defect of detection is noted or suspected, the lens should be returned to Hanita Lenses.

For more quality ophthalmic products, Implant & Folding instructions, please log onto our web site at www.hanitалenses.com

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