

MONOPLUS · ENHANCED INTERMEDIATE

---

# VISION EXTENDED




Extend™

World's First  
Bessel-Optics Monofocal Plus IOL

HANITA  
Lenses


# The Extend™ Monofocal Plus IOL




Bessel beams, characterized by their unique non-diffractive propagation, "needle-shaped beam"<sup>[1]</sup>, and self-healing capabilities, are increasingly applied in ophthalmology for high-resolution imaging and enhanced fixation targets<sup>[2]</sup>.

## Clinical Advantages from Bessel Optics

The Extend™ IOL, powered by Bessel optics, delivers an extended depth of focus, rather than a single focal point, to cover both distance and intermediate vision without the complexity of multifocal designs.




The self-healing mechanism<sup>[3]</sup> of Bessel beams enables them to restore their structure when partially obstructed, as unblocked wavefronts re-interfere to maintain stable imaging. This suggests the IOL might tolerate optical path obstructions like low-level posterior capsule opacification (PCO) or floaters.



1. Zhao, J., Winetraub, Y., Du, L., VAN Vleek, A., Ichimura, K., Huang, C., Aasl, S. Z., Sarin, K. Y., & de la Zerda, A. (2022, August 20). Flexible method for generating needle-shaped beams and its application in optical coherence tomography. *Optica*. <https://pmc.ncbi.nlm.nih.gov/articles/PMC10243785/>
2. Suchand Sandeep, C. S., Khairyanto, A., Aung, T., & Vadakke Matham, M. (2023). Bessel beams in ophthalmology: A Review. *Micromachines*, 14(9), 1672. <https://doi.org/10.3390/mi14091672>
3. Aiello, A. (2014, December). Wave-optics description of self-healing mechanism in Bessel beams. ResearchGate. [https://www.researchgate.net/publication/266971869\\_Wave-optics\\_description\\_of\\_self-healing\\_mechanism\\_in\\_Bessel\\_beams](https://www.researchgate.net/publication/266971869_Wave-optics_description_of_self-healing_mechanism_in_Bessel_beams)

# Extend Your Patients' Vivid Moments With Clarity

Distance vision comparable to a standard monofocal IOL, along with functional intermediate vision.



Extend™: Bianchi et al., Highlights of Ophthalmology 2023 Volume 51(2)

## Full Range Binocular Vision


By using a mini-monovision approach targeting small myopia in the non-dominant eye, a full range of functional vision can potentially be achieved.

One eye: **emmetropia**


The other eye: **mild myopia**




# Extend™: Bessel-Powered Monofocal Reinvented




When tested with USAF target images, monofocal IOLs are sharp only at their preset focal point (e.g., distance "far-point" row), with marked blurring at off-focus states (-1.0D to -1.5D). In contrast, Extend™ IOL, using Bessel-powered extended focus technology, maintains discernible target clarity across all tested refractive states, enabling continuous distance-to-intermediate vision and outperforming monofocal designs.







Polychromatic PSF images taken at optimal saturation with a pupil diameter of 4.5 mm. Both IOLs showed comparable light distribution with a speckle-like pattern outside the PSF core.

Simulated VA demonstrated that at far focus, the Extend™ IOL achieves a high-acuity logMAR of -0.09, with stable vision maintained across a broad defocus range, enabling clear sight from far to intermediate distances.




The Extend™ IOL delivers strong MTF performance: high detail retention at 50 lp/mm (0.32) with a 4.5-mm pupil for low-light clarity, and reliable contrast in bright 3-mm pupil scenarios.<sup>[1]</sup>



1. Łabuz, G., Yan, W., & Munro, D. J. (2021). The optical-quality assessment of the Hanita mono-EDoF intraocular lens versus a competitor model [Unpublished report]. David J Apple International Laboratory for Ocular Pathology, University of Heidelberg.



By addressing three key aspects—light-spot focusing, vision stability, and detail recognition in both day and night conditions, data show that the Hanita Lenses Extend™ IOL provides sharp distance vision, covers intermediate viewing ranges, and maintains clarity in low-light environments. This makes it particularly beneficial for patients who require high-quality postoperative vision, such as for night driving.



# Hanita Lenses Extend™ SL

## OPTIC CHARACTERISTICS

Diopter range (SE)	+5.0 to +34.0 D
Diopter increment	0.5 D (+5.0 to +30.0 D) 1.0 D (+30.0 to +34.0 D)
Optical design	Modified high-order aspheric for extension of the depth of focus

## GEOMETRIES

Optic diameter	6.0 mm
Total diameter	13.0 mm
Haptic configuration	C-loop
Edge design	Square edge
Haptic angulation	5°

## MATERIAL SPECIFICATIONS

Material	Hydrophobic acrylic with bonded UV absorber and violet light filter
Refractive index	1.48 (@ 35 °C)
Abbe number	49
Filtration	UV and violet light

## OPTICAL BIOMETRY

SRK/T: A-constant	119
Holladay II: ACD	5.549


## CONTACT ULTRASOUND BIOMETRY

SRK/T: A-constant	118.5
Holladay II: ACD	5.257

## RECOMMENDED DELIVERY SYSTEM



Preloaded  
Injector: Accuject 2.2



# Hanita Lenses Extend™ SL Toric

## OPTIC CHARACTERISTICS

Diopter range (SE)	+5.0 to +34.0 D
Diopter increment	0.5 D (+5.0 to +30.0 D) 1.0 D (+30.0 to +34.0 D)
Cylinder range	1.0, 1.5, 2.25, 3.0, 3.75, 4.5 D
Optical design	Modified high-order aspheric for extension of the depth of focus

## GEOMETRIES

Optic diameter	6.0 mm
Total diameter	13.0 mm
Haptic configuration	C-loop
Edge design	Square edge
Haptic angulation	5°

## MATERIAL SPECIFICATIONS

Material	Hydrophobic acrylic with bonded UV absorber and violet light filter
Refractive index	1.48 (@ 35 °C)
Abbe number	49
Filtration	UV and violet light

## OPTICAL BIOMETRY


SRK/T: A-constant	119
Holladay II: ACD	5.549

## CONTACT ULTRASOUND BIOMETRY

SRK/T: A-constant	118.5
Holladay II: ACD	5.257

## CALCULATOR

## RECOMMENDED DELIVERY SYSTEM



Preloaded  
Injector: Accuject 2.2

# Extend™ IOL by Hanita Lenses

Pioneering a new era  
of monofocal IOL technology.

- Distance vision comparable to standard monofocal IOLs
- Functional intermediate vision
- Excels in low light
- More tolerant of visual disturbances
- Overall performance comparable to the industry benchmark



For healthcare professionals only.

Please read the Instruction for Use for important safety information and consult our specialists if you have any questions.



**Extend™**

Extend™ is a trademark of Hanita Lenses  
© 2025 Hanita Lenses | [www.hanitalenses.com](http://www.hanitalenses.com)

The information contained herein is subject to change without notice.

Hanita Lenses shall not be liable for technical or editorial errors or omissions contained herein.

GRA-0000770 Rev. 3.0