

MONOFOCAL PLUS • ENHANCED INTERMEDIATE

VISION EXTENDED




Extend™

World's First
Bessel-Optics Monofocal Plus IOL

HANITA
Lenses


The Extend™ Monofocal Plus IOL




Bessel beams, characterized by their unique non-diffractive propagation, "needle-shaped beam" profile^[1], and self-healing capabilities, are increasingly applied in ophthalmology for high-resolution imaging and enhanced fixation targets^[2].

Clinical Advantages from Bessel Optics

The Extend™ IOL, powered by Bessel optics, delivers an extended depth of focus rather than a single focal point, to cover both distance and intermediate vision without the complexity of multifocal designs.




The self-healing mechanism^[3] of Bessel beams enables them to restore their structure when partially obstructed, as unblocked wavefronts re-interfere to maintain stable imaging. This suggests the IOL might tolerate optical path obstructions like low-level posterior capsule opacification (PCO) or floaters.



1. Zhao, J., Winetraub, Y., Du, L., Van Vleck, A., Ichimura, K., Huang, C., Aasi, S. Z., Sarin, K. Y., & de la Zerna, A. (2022, August 20). Flexible method for generating needle-shaped beams and its application in optical coherence tomography. *Optica*. <https://pmc.ncbi.nlm.nih.gov/articles/PMC10243785/>
2. Suchand Sandeep, C. S., Khairyanto, A., Aung, T., & Vadakke Matham, M. (2023). Bessel beams in ophthalmology: A review. *Micromachines*, 14(9), 1672. <https://doi.org/10.3390/mi14091672>
3. Aiello, A. (2014, December). Wave-optics description of self-healing mechanism in Bessel beams. ResearchGate. https://www.researchgate.net/publication/266971869_Wave-optics_description_of_self-healing_mechanism_in_Bessel_beams

Extend Your Patients' Vivid Moments With Clarity

Distance vision comparable to a standard monofocal IOL, along with functional intermediate vision.



Extend™: Bianchi et al., Highlights of Ophthalmology 2023;51(2).

Full Range Binocular Vision


By using a mini-monovision approach targeting low myopia in the non-dominant eye, a full range of functional vision can potentially be achieved.

One eye: **emmetropia**


The other eye: **mild myopia**




Extend[™]: Bessel-Powered Monofocal Reinvented




When tested with USAF target images, monofocal IOLs are sharp only at their preset focal point (e.g., distance "Far-point" row), with marked blurring at off-focus states (-1.0 D to -1.5 D). In contrast, the Extend[™] IOL, using Bessel-powered extended-focus technology, maintains discernible target clarity across all tested refractive states, enabling continuous distance-to-intermediate vision and outperforming monofocal designs.







Polychromatic PSF images taken at optimal saturation with a pupil diameter of 4.5 mm. Both IOLs show comparable light distribution with a speckle-like pattern outside the PSF core.

Simulated VA demonstrated that at far focus, the Extend™ IOL achieves a high-acuity logMAR of -0.09, with stable vision maintained across a broad defocus range, enabling clear sight from far to intermediate distances.




The Extend™ IOL delivers strong MTF performance: high detail retention at 50 lp/mm (0.32) with a 4.5-mm pupil for low-light clarity, and reliable contrast in bright-light conditions with a 3-mm pupil^[1].



1. Kremser, F., Yan, W., Khoramnia, R., et al. Monofocal intraocular lens based on the Bessel principle for improved intermediate vision: a comparative assessment. BMC Ophthalmol 26, 336 (2026). <https://doi.org/10.1186/s12886-026-04973-9>



By addressing three key aspects—light-spot focusing, vision stability, and detail recognition in both day and night conditions—data show that the Hanita Lenses Extend™ IOL provides sharp distance vision, covers intermediate viewing ranges, and maintains clarity in low-light environments. This makes it particularly beneficial for patients who require high-quality postoperative vision, such as night driving.



Hanita Lenses Extend™ SL

OPTIC CHARACTERISTICS

Diopter range (SE)	+5.0 to +30.0 D
Diopter increment	0.5 D
Optical design	Modified high-order aspheric for extension of the depth of focus

GEOMETRIES

Optic diameter	6.0 mm
Total diameter	13.0 mm
Haptic configuration	C-loop
Edge design	Square edge
Haptic angulation	5°

MATERIAL SPECIFICATIONS

Material	Hydrophobic acrylic with bonded UV absorber and violet light filter
Refractive index	1.48 (@ 35 °C)
Abbe number	49
Filtration	UV and violet light

OPTICAL BIOMETRY

SRK/T: A-constant	119
Holladay II: ACD	5.549


CONTACT ULTRASOUND BIOMETRY

SRK/T: A-constant	118.5
Holladay II: ACD	5.257

RECOMMENDED DELIVERY SYSTEM



Preloaded
injector: Accuject 2.2



Hanita Lenses Extend™ SL Toric

OPTIC CHARACTERISTICS

Diopter range (SE)	+5.0 to +27.5 D
Diopter increment	0.5 D
Cylinder range	1.0, 1.5, 2.25, 3.0 D
Optical design	Modified high-order aspheric for extension of the depth of focus

GEOMETRIES

Optic diameter	6.0 mm
Total diameter	13.0 mm
Haptic configuration	C-loop
Edge design	Square edge
Haptic angulation	5°

MATERIAL SPECIFICATIONS

Material	Hydrophobic acrylic with bonded UV absorber and violet light filter
Refractive index	1.48 (@ 35 °C)
Abbe number	49
Filtration	UV and violet light

OPTICAL BIOMETRY

SRK/T: A-constant	119
Holladay II: ACD	5.549

CONTACT ULTRASOUND BIOMETRY

SRK/T: A-constant	118.5
Holladay II: ACD	5.257

CALCULATOR

RECOMMENDED DELIVERY SYSTEM

K Notation: ☐ Dpt ☐ Wave

OD ☐ OS

43 @ 170
Flat K (dpt) Axis (°)

45 @ 80
Steep K (dpt) Axis (°)


1.3375 1.332 1.3315
n-value

23 3.2
Axial Length (mm) ACD (optional) (mm)

0.2 @ 90
SIA (dpt)

☒ Include posterior corneal astigmatism

Calculate



Preloaded
injector: Accuject 2.2

Extend™ IOL by Hanita Lenses

Pioneering a new era
of monofocal IOL technology.

- Distance vision comparable to standard monofocal IOLs
- Functional intermediate vision
- Excels in low light
- More tolerant of visual disturbances
- Overall performance comparable to the industry benchmark



For healthcare professionals only.

Please read the Instructions for Use for important safety information and consult our specialists if you have any questions.



Extend™

Extend™ is a trademark of Hanita Lenses.
© 2026 Hanita Lenses | www.hanitalenses.com

The information contained herein is subject to change without notice.

Hanita Lenses shall not be liable for technical or editorial errors or omissions contained herein.

GRA-0000770 Rev. 4.0