



# Monofocal IOLs

## Aspherical HP IOLs

\* SeeLens HP



\* BunnyLens HP



\* Available in preloaded configuration

	SeeLens HP *	BunnyLens HP *
<b>Platform</b>	C-Loop	4 - Loop
<b>Overall Diameter</b>	13.0mm	11 mm (>10D) 11.5mm (≤10D)
<b>Placement</b>	Capsular Bag	
<b>Optic diameter</b>	6.0 mm	
<b>Power range</b>	+5.0 to +10.0 (1D increments) +10.0 to +30.0 (0.5D increments) +30.0 to +35.0 (1D increments)	
<b>Optic design</b>	Aspherical	
<b>360° Continuous Square Edge</b>	Yes	
<b>Haptic angulation</b>	5°	
<b>Material</b>	Hydrophobic acrylic with bonded UV absorber and violet light filter	
<b>Refractive Index</b>	1.48 (@ 35°C)	
<b>A-constant (SKR/T) for Optical or Immersion US biometry</b>	119.5	118.5
<b>A-constant (SKR/T) for Contact US biometry</b>	119.16	118.16
<b>Sterilization</b>	ETO	

# Providing Excellent Vision Quality For Your Patient

The Aspherical HP IOLs were designed by the most advanced optical tools to provide:

- reduced Spherical aberration
- functional vision
- vision in night conditions

The Aspherical Hydrophobic lens from Hanita Lenses, provides the patient with an excellent vision quality at day and night conditions, by using state-of-the-art aberration free Aspherical optical design.

## Stability and Centration

