


## Technical data sheet

Trade name: **Softject 2.4 – 1P**

Designation: Intra ocular lens Injector

Description:

- a sub assembly: injector with a body, a pusher, a metal spring and a plunger.
- a cartridge Medical Ref. LC 604251
- An accessory: Reamer



Raw material

(MABS) Methyl Methacrylate Acrylonitrile Butadiene Styrene, (POM) Polyoxymethylene, (SEBS) Styrene Ethylene Butadiene Styrene, (PP) Polypropylene, (inox steel AISI 302).  
Latex free

Sterilization process: ETO

|                     |  |
|---------------------|--|
| Reference<br>HANITA | XPH VG 2.2   |
| Packaging           | 2 boxes of 50 devices in blister grouped in a carton |

T1845 03/2013